



#### Finishes Available

- |   |                   |   |        |
|---|-------------------|---|--------|
| 1 | Natural Oiled     | 4 | Olive  |
| 2 | Natural Lacquered | 5 | Smoked |
| 3 | Lime              | 6 | Milan  |



Individual options available



## Theramatic - solid natural oak

Handcrafted in beautiful solid oak, these stylish beds bear no resemblance to a hospital type bed in looks, but have all of the assistive function that you need in helping to maintain independence at home.

Your bed is specified to your exact requirements meaning that you can have any combination of functions and sizes together in a bed that looks like part of your home.

Why not choose a beautiful solid oak bed enclosure to surround the metal mechanism that powers a Theraposture bed? There is a choice of wood finishes and colours:

Also you can choose between a slatted or panelled bed surround with high or low foot ends.

We also provide an extensive range of matching solid oak bedroom furniture such as bedside cabinets, chests of drawers and wardrobes. For further examples from this impressive collection, see our accessories page.

### Oak beds for individuals

Theraposture's stylish solid oak beds look at home in your home.

They also help overcome everyday difficulties through their functionality and provide you with the opportunity to remain independent in your own home without the sacrifice of intrusive adaptation.

Beds for individuals can be any size and can have any function to help with everyday positioning and transfer difficulties. (See our technical data for function and size options).

### Theramatic mechanisms - Profiling only

The profiling only action will allow you to independently raise and lower both the back rest and the leg section by simply pushing a button on the hand controller. This mechanism is often the choice of partners in a double combination bed as it allows both individuals to sit up independently together. The metal construction makes this mechanism more durable for long term use and the only choice if grab handles or side rails are needed.

### Variable height and profiling

In addition to the profiling only mechanism the variable height function is of significant benefit to less able people in a number of ways. If a carer or family member assists a person in bed, the variable height function will smoothly and safely raise the occupant to a safe working height, thereby reducing the risk of a back care injury.

Additionally, this function will help a person to achieve a standing transfer by vertically raising them to standing from sitting in a controlled and safe movement.

These strong steel mechanisms are exclusive to Theraposture. Many of our beds are used by clients with long term conditions and we developed these mechanisms, in consultation with our clients to last a very long time. The steel makes the mechanism very strong and the actuators and other electrics used are of the highest grade. These beds are made for the higher usage that many of our clients demand.



Variable Height and Profiling mechanism



Profiling only mechanism



Easy to use hand control with contrasting buttons

**Freephone: 0800 834654**

Email: [info@theraposture.co.uk](mailto:info@theraposture.co.uk) [www.theraposture.co.uk](http://www.theraposture.co.uk)

Theraposture Limited, Kingdom Avenue, Northacre Industrial Park, Westbury, Wiltshire. BA13 4WE  
Open: Monday – Friday, 8:30am–5:30pm. Answerphone messages can be left at all other times

